

Administer Mode



"Check Before You Inject"

An additional safety check, **Administer Mode** allows a Safe Label System (SLS) user to **"check before you inject"** to confirm a medication prior to administration to a patient. When a medication is being prepared, SLS audibly states and visually displays the drug name and concentration and immediately prints a ready-to-apply TJC and ASA compliant color label. Once that medication is prepared and labeled, the barcode on the SLS label enables a "triple check" using Administer Mode. When the clinician is ready to administer the medication, they simply scan the labeled medication and the SLS Point of Care Station (PCS) presents another audible and visual confirmation of the drug name and concentration, as well as displays the time remaining until the medication will expire. Administer Mode greatly reduces the risk of syringe swaps, one of the three most common medication errors occurring in the operating room. With Administer Mode enabled, the Safe Label System fully meets the APSF recommendations to scan drugs before preparation and/or administration.

SLS PCS showing
Administer Mode



KEY FEATURES

Advantages of Administer Mode

- Provides a "triple check" of medication prior to administration to the patient
- Meets APSF recommendations to scan drugs before preparation and/or administration
- Visually and audibly confirms the drug, concentration and time until expiration
- Greatly reduces the risk of syringe swaps
- Optional hand scanner enables the preparation and administration of items unable to be put under the SLS scanner, whether due to proximity or size

Call Codonics today at **800.444.1198** or visit **www.codonics.com** for more information.

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents: www.codonics.com/ip/patents Printed in the U.S.A. Copyright © 2016 Codonics Inc. 9/2016



Sponsoring Donor



Corporate Member



We bring the future into focus

17991 Englewood Drive
Middleburg Heights, OH 44130 USA
+1.440.243.1198
+1.440.243.1334 Fax
info@codonics.com
www.codonics.com